



VoltaFeet:

Verification and Validation

Senior Design Group 2

Ashton Naumann, Logan Groneck, Matt Heiken

In collaboration with Dr. Dan Moran

Presented by Matt Heiken

Project Overview



Goal:


To build a **simple, long lasting**, treatment for **eliminating** toe/foot **microbes** as well as shoe and foot **odor** for a diverse population.

Project Overview



Goal:


To build a **simple, long lasting**, treatment for **eliminating** toe/foot **microbes** as well as shoe and foot **odor** for a diverse population.



Users with:


- Bromodosis (stinky feet)
- Athlete's Foot
- Foot Wounds
- Diabetic Ulcers

Project Overview



A sock that incorporates a mesh of **silver wire** and a **3V lithium cell battery** with a **constant current source** to conduct microamps of **DC current** for **continuous and controlled release** of antimicrobial **silver ions** onto the surface of the foot.

Project Overview



A sock that incorporates a **conductive silver fabric** and a **3V lithium cell battery** with a **constant current source** to conduct microamps of **DC current** for **continuous and controlled release** of antimicrobial **silver ions** onto the surface of the foot.

Updated Specifications



Category	Old Spec	New Spec
Diverse Sizing	Reliable for feet sized 4 through 12 men's and women's.	Manufactured in three different sizings (small, medium, and large) that encompasses feet sized 4 through 14 men's and women's
User Safety	Must not result in shocks or punctures after repeated loading during locomotion in any practical environment; If an electrical current is used, it must be less than 1 mA.	Must not result in shocks or punctures after repeated loading during locomotion in any practical environment; If an electrical current is used, it must be less than 5 mA.

New Specification



Category	Specification
Antimicrobial Penetration	Silver ion depth penetration comparable to other antimicrobial silver treatments.

New Specification



Category

Specification

**Antimicrobial
Penetration**

Silver ion depth penetration comparable to other antimicrobial silver treatments.



No changes to need statement, project scope, design schedule, and team responsibilities



1





Verification Plan

Was the final product built correctly?

Important Considerations





Important Considerations




Grid size optimization: Finding the largest grid spacing that still provides antimicrobial effects to the entire area





Important Considerations



Grid size optimization: Finding the largest grid spacing that still provides antimicrobial effects to the entire area





Important Considerations





Silver penetration deeper than **0.1 mm**

Important Considerations



Current optimization: Balancing antimicrobial efficiency with battery life

Important Considerations




Eliminates 99.9% of microbes/fungus in **12 hours**.

3D Cell Culture




- 3D soft agar E. coli culture
- Independent variables
 - Current amplitude
 - Surface Area
 - Type of Fabric/Mesh
- Measure **depth and coverage** over time



2D Cell Culture




- Testing the efficacy of silver on **multiple different microbial and fungal strains**



2D Cell Culture




- Testing the efficacy of silver on **multiple different microbial and fungal strains**
 - **However...** difficulties





Validation Plan

Did we build the right product?



Important Considerations



**Sterile
Foot
Surface**



Comfort




**Functional
Duration**



Weight

Testing Foot Sterility




Bilateral Control
Experiment

Testing Foot Sterility



Bilateral Control
Experiment



At different points in time:

- Foot swabs spread on nutrient rich agar plate and incubated from different parts of the foot

Independent Variable: Sock type

Dependent Variable: Colony count

Validating Quality



Face-to-Face
Survey

Comparison of commercial prototype to normal sock

- Likert Scale (Ratings ranging from 1-5)
- Tested Specifications
 - Comfort
 - Appearance



Testing Longevity

> month

Limiting Factor: Battery vs. Silver Supply

Battery

Wearing the device on a regular basis and **frequently measuring voltage** with a voltmeter.


vs.

Silver Supply

Impedance curves: Measuring the impedance of the conductive fabric over time to determine the loss of silver

Testing Longevity

> month



Silver Supply

Impedance curves: Measuring the impedance of the conductive fabric over time to determine the loss of silver

3



FDA Approval

How will we get our device approved?

VoltaFeet will rely on 510(k) premarket notification

“ A 510(K) is a premarket submission made to FDA to demonstrate that the device to be marketed is **at least as safe and effective**, that is, substantially equivalent, **to a legally marketed device** that is not subject to premarket approval. ”

-FDA.gov

VoltaFeet will rely on 510(k) premarket notification

“ A 510(K) is a premarket submission made to FDA to demonstrate that the device to be marketed is **at least as safe and effective**, that is, substantially equivalent, to a **legally marketed device** that is not subject to premarket approval. ”

-FDA.gov

Legally Marketed Devices

- Acticoat Flex 7 Dressing (K083113) - **Class II**
 - Uses **silver coated nylon** and **silver ions** for **antimicrobial activity on the surface of the skin**
- Transcutaneous Electrical Nerve Stimulation (TENS) Device - **Class II**
 - utilizes **0-50 mA DC current on skin**





Current Status of Project




Approval for cell culture



Current Status of Project



Approval for cell culture



Proof of Concept
Testing




Current Status of Project



Approval for cell culture



Proof of Concept
Testing




Prototyping

Current Status of Pro



Approval for cell culture




Prototyping



Current Status of P



Approval for cell culture



Prototyping





Current Status of Project



Approval for cell culture



Prototyping

→ How do we do this better?



Thank you!

Questions?